

# E Series

Handheld UV Lamp

ENF-240C + ENF-260C + ENF-280C



Forensics



Lab &  
Life Sciences



Mineral  
Inspection

**ERGONOMIC CARRYING HANDLE**  
Comfortable grip while operating the lamp.

**TOGGLE SWITCH**  
Built in toggle switch for easy On/Off operation



**ANODIZED ALUMINUM HOUSING**  
Lightweight and rugged for the toughest work environments.

**LONG LIFE FILTER GLASS**  
Our unique filter glass ages at a slower rate than conventional filters and resists solarization extending the lifetime of the lamp.





**ENF-240C E SERIES HANDHELD UV LAMP**

**ENF-260C E SERIES HANDHELD UV LAMP**



**ENF-280C E SERIES HANDHELD UV LAMP**

	ENF-240C	ENF-260C	ENF-280C
<b>Light Source:</b>	4 Watt 365nm Tube / 4 Watt 254nm Tube with filter assembly	4 Watt 365nm Tube / 6 Watt 254nm Tube with filter assembly	8 Watt 365nm Tube / 8 Watt 254nm Tube with filter assembly
<b>Lamp Style:</b>	Handheld Lamp	Handheld Lamp	Handheld Lamp
<b>Length:</b>	8 in (20.32 cm)	10.25 in (26.04 cm)	13.25 in (33.65 cm)
<b>Width:</b>	3.25 in (8.26 cm)	3.25 in (8.26cm)	3.25 in (8.26cm)
<b>Height:</b>	2.5 in (6.35 cm)	2.5 in (6.35 cm)	2.5 in (6.35cm)
<b>Weight:</b>	2.82 lb (1.28 kg)	3.9 lb (1.77 kg)	1.6 lb (.72 kg)
<b>Plug Options:</b>	USA, Australian (/FA), United Kingdom (/FBE), China (/FC), European (/FE), Japan(/J), 12 Volt (/12H)		