

# MINI LAMP

Battery Operated UV Lamp

UV-4B



Forensic Investigation



Lab & Life Sciences



Rodent Detection

## POWER SWITCH

Textured On/Off switch for easy one handed operation.

## HOUSING

Made from high quality ultra light ABS plastic and nylon for fatigue free use.

## 4 WATT SELF FILTERED UV TUBE

Lamp reaches full intensity immediately with no warm up period required.

## POWER

Powered by four (4) AA (1.5 Volt) alkaline batteries (not included). Lamp features a sliding cover for easy battery replacement.

## NYLON LANYARD

Includes nylon lanyard to secure handheld UV lamp to users hand while in use.





<b>Light Source:</b>	(1) 4 Watt UV-A Tube
<b>Lamp Style:</b>	Handheld Lamp
<b>Dimensions:</b>	6 ¼ x 2 x 1in (15.9 x 5.1 x 2.5 cm)
<b>Weight:</b>	7 oz (198 g) (with Battery)
<b>Power Requirements:</b>	(4) AA (1.5 Volt) Alkaline Batteries
<b>Run Time:</b>	4 to 6 hours (continuous)

**Blacklight Lamp Specialty Applications**

- Counterfeit Bill Detection
- Detecting Pet Urine
- Rodent Detection
- Pest Control
- Hidden Stains
- Inspecting Stamps
- Mineral and Gem Inspection
- Fine Art and Antique Inspection: Detecting porcelain & pottery repairs, testing for reproduction glass, examining artwork, dating textiles & fabrics, authenticating cast iron
- Document Verification: Confidential documents, checks, tickets, passports
- Watermarks, stains, imperfections