

B-260

Battery-Operated UV Lamp

SPECTRO-UV®

Rechargeable Handheld UV Blacklight for Professional Applications.

The B-260 rechargeable handheld UV black-light lamp offers a high-performance solution at an economical price. This portable UV lamp features a dual Ultraviolet (UV-A) 365 nm tube design, a contoured carrying handle and a 6-Volt rechargeable battery system. Constructed from durable, high-impact plastic housing, the B-260 ensures years of reliable use in the toughest work environments.



Key Features

- **Dual-Tube Performance:** Dual 6-watt UV-A (365nm) tubes deliver optimal fluorescence detection
- **Extended Battery Operation:** 6-volt, 4 amp-hour rechargeable battery provides extended field use
- **Professional Construction:** Durable, high-impact plastic housing withstands demanding work environments
- **User-Friendly Design:** Built-in charging indicator (red LED) shows charging status
- **Economical price point** with professional performance
- **Compact, portable design** for easy transport
- **Built-in charger** eliminates additional accessories

Part Number	Light Source	Dimensions	Weight	Power	Charger	Charging Time
B-260	6-Watt UV-A tube	10.25" x 3.15" x 5.63" (26 x 8 x 14.3 cm)	Lightweight portable design	6-volt, 4 amp-hour rechargeable battery	Built-in AC charger (120 V or 230 V)	Approx. 24 hours (initial) 16-24 hours (subsequent)

SPECTRO-UV®

www.Spectro-UV.com
4 Dubon Ct., Farmingdale, NY 11735
866-230-7305

LAB
12/25 AB260.1
PRINTED IN USA
©2025. Spectro-UV. All rights reserved.