

SPECTROLINE®

CM Series

Fluorescence Analysis Cabinet

CM-10A + CM-26A



Forensics



Lab &
Life Sciences



Mineral
Inspection

CONTOURED EYEPiece

Flexible, contoured eyepiece with a built in UV absorbing window for safe sample viewing.

UV LAMP (MODEL CM-26A ONLY)

Outfitted with removable UV light source producing a long wave (365nm) intensity of 300 $\mu\text{W}/\text{cm}^2$ and a shortwave (254nm) intensity of 310 $\mu\text{W}/\text{cm}^2$ at the cabinet floor.



ALUMINUM HOUSING

FELT CURTAINS

Felt curtains permit easy access to cabinet interior while shielding samples from external light

REMOVABLE BOTTOM PANEL

Removable bottom panel to place cabinet over other devices and equipment.



CM-10A

Length:	12 in (30 cm)
Width:	9 in (22.8 cm)
Height:	6.5 in (16.5 cm)
Weight:	7 lbs(3.2kg) with ENF-260C UV Lamp



CM-26A

Length:	12 in (30 cm)
Width:	9 in (22.8 cm)
Height:	6.5 in (16.5 cm)
Weight:	7 lbs(3.2kg) with ENF-260C UV Lamp
Light Source:	ENF-260C
Lamp Style:	Handheld UV Lamp



SPECTRO-UV

www.Spectro-UV.com

4 Dubon Ct., Farmingdale, NY 11735
866-230-7305

DISTRIBUTED BY: