

E SERIES HANDHELD LAMP MATRIX

Part Number	Light Source	Lamp Style	Length	Width	Height	Weight	Applications
E Series							
ENF-240C	4 Watt 365nm Tube / 4 Watt 254nm Tube with filter assembly	Handheld Lamp	8 in (20.32 cm)	3.25 in (8.26 cm)	2.5 in (6.35 cm)	2.82 lb (1.28 kg)	Lab & Life Sciences, Agriculture/Rodent Contamination Detection, Conservation, Gem & Mineral Inspection
ENF-260C	6 Watt 365nm Tube / 6 Watt 254nm Tube with filter assembly	Handheld Lamp	10.25 in (26.04 cm)	3.25 in (8.26 cm)	2.5 in (6.35 cm)	3.9 lb (1.77 kg)	Lab & Life Sciences, Agriculture/Rodent Contamination Detection, Conservation, Gem & Mineral Inspection
ENF-280C	8 Watt 365nm Tube / 8 Watt 254nm Tube with filter assembly	Handheld Lamp	13.25 in (33.65 cm)	3.25 in (8.26 cm)	2.5 in (6.35 cm)	1.6 lb (.72 kg)	Lab & Life Sciences, Agriculture/Rodent Contamination Detection, Conservation, Gem & Mineral Inspection
EA-140	UV-A Tube, 4 Watts	Handheld Lamp	7 in (17.8 cm)	3.25 in (8.26 cm)	2.5 in (6.35 cm)	2.45 lb (1.1 kg)	Lab & Life Sciences, Agriculture/Rodent Contamination Detection, Conservation, Gem & Mineral Inspection
EA-160	UV-A Tube, 6 Watts	Handheld Lamp	10.25 in (26.04 cm)	3.25 in (8.26 cm)	2.5 in (6.35 cm)	2.7 lb (1.22 kg)	Lab & Life Sciences, Agriculture/Rodent Contamination Detection, Conservation, Gem & Mineral Inspection

Key Selling Points:

- Choice of 4-, 6- and 8- watt models
- Short wave UV-C lamps are ideal for RUVIS
- Integrally filtered tubes or external filter assemblies
- Durable, silver-anodized aluminum housings

