

X-SERIES BENCH LAMP MATRIX

Part Number	Light Source	Lamp Style	Power	Wavelength	No of Tubes	Length	Width	Height	Weight	Applications
X-Series										
X-15G	UV-C Tube	Bench Lamp	15 Watts	UV-C (254 nm)	1	18.75 in (47.6 cm)	4.5 in (11.4 cm)	3.5 in (8.9 cm)	6.11 lb (2.77 kg)	Lab & Life Sciences, Conservation, NDT
XX-15A	UV-A Tube	Bench Lamp	15 Watts	UV-A (365 nm)	2	18.75 in (47.6 cm)	4.5 in (11.4 cm)	3.5 in (8.9 cm)	6 lb (2.72 kg)	Lab & Life Sciences, Conservation, NDT
XX-15BF	UV-B/UV-C Tube	Bench Lamp	15 Watts	UV-B (312 nm) / UV-C (254 nm)	1/1	18.75 in (47.6 cm)	4.5 in (11.4 cm)	3.5 in (8.9 cm)	7.32 lb (3.32 kg)	Lab & Life Sciences, Conservation, NDT
XX-15G	UV-C Tube	Bench Lamp	15 Watts	UV-C (254 nm)	2	18.75 in (47.6 cm)	4.5 in (11.4 cm)	3.5 in (8.9 cm)	6.11 lb (2.77 kg)	Lab & Life Sciences, Conservation, NDT
XX-40	UV-A Tube	Bench Lamp	40 Watts	UV-A (365 nm)	2	48 in (121.9 cm)	6 in (15.1 cm)	4 in (10.2 cm)	17 lb (7.71 kg)	Lab & Life Sciences, Conservation, NDT

Key Selling Points:

- Choice of 15- or 40-watt models
- Available in all UV wavelengths (LW, MW and SW)
- Filtered UV-A and unfiltered UV tubes
- Offered with our without filter assemblies

