

CROSSLINKER MATRIX

Part Number	Overall Housing Dimensions (W x H x D)	Effective Inner Chamber Dimensions (W x H x D)	Door Dimensions (W x H)	Viewing Window Dimensions (W x H)	Net Weight	Power	Wave-length	Tubes	Applications
XL-1000	19.5" x 10.5" x 9" (49.5 x 26.7 x 22.9 cm)	13.5" x 7" x 7.5" (34.3 x 17.8 x 19.1 cm)	12" x 6.25" (30.5 x 15.9 cm)	6" x 2.5" (15.2 x 6.4 cm)	17.5 lb (7.9 kg)	120V-60Hz-2A, 230V-50Hz-2A, 100V-50/60Hz-2A	254nm	(5) 8W SW, BLE-8T254	Lab & Life Sciences, Forensics
XL-1500	24" x 10.5" x 14" (61 x 26.7 x 35.6 cm)	18.25" x 6.25" x 12.5" (46.4 x 15.9 x 31.8 cm)	16" x 6.25" (40.6 x 15.9 cm)	6" x 2.5" (15.2 x 6.4 cm)	26 lb (11.8 kg)	120V-60Hz-3A, 230V-50Hz-3A, 100V-50/60Hz-3A	254nm	(6) 15W SW, BLE-1T155	Lab & Life Sciences, Forensics

Key Selling Points:

- Rugged, anodized aluminum housing
- Long lasting anodized aluminum chamber prevents rusting
- Available in SW, MW and LW wavelengths
- Available in 2 sizes: Standard and Large

