

E SERIES

Handheld UV Lamps

SPECTRO-UV[®]

Professional-grade handheld UV lamps engineered for laboratory and field applications.

The **E Series** delivers dual-wavelength capability (254 nm/365 nm) in three power configurations, featuring rugged anodized aluminum construction and our exclusive long-life filter technology.



Key Features

- **Dual Wavelength Options:** 254nm and 365nm tubes with precision filter assemblies
- **Three Power Levels:** 4W, 6W and 8W configurations to match your application needs
- **Ergonomic Design:** Comfortable carrying handle reduces operator fatigue
- **Rugged Construction:** Lightweight anodized housing withstands demanding environments
- **Long-Life Filter Glass:** Proprietary filter technology resists solarization and extends lamp lifetime
- **Simple Operation:** Built-in toggle switch for easy on/off control
- **Global Compatibility:** Multiple plug options available for world-wide use



ENF-240C



ENF-260C



ENF-280C



EA-140



EA-160

Part Number	Light Source	Length	Width	Height	Weight	Plug Options
ENF-240C	4 Watt 365nm Tube / 4 Watt 254nm Tube with filter assembly	8 in (20.32 cm)	3.25 in (8.26 cm)	2.5 in (6.35 cm)	2.82 lb (1.28 kg)	USA, Australian, United Kingdom, China, Europe, Japan, 12V
ENF-260C	6 Watt 365nm Tube / 6 Watt 254nm Tube with filter assembly	10.25 in (26.04 cm)	3.25 in (8.26 cm)	2.5 in (6.35 cm)	3.9 lb (1.77 kg)	USA, Australian, United Kingdom, China, Europe, Japan, 12V
ENF-280C	8 Watt 365nm Tube / 8 Watt 254nm Tube with filter assembly	13.25 in (33.65 cm)	3.25 in (8.26 cm)	2.5 in (6.35 cm)	1.6 lb (.72 kg)	USA, Australian, United Kingdom, China, Europe, Japan, 12V
EA-140	UV-A Tube, 4 Watts	7 in (17.8 cm)	3.25 in (8.26 cm)	2.5 in (6.35 cm)	2.45 lb (1.1 kg)	USA, Australian, United Kingdom, China, Europe, Japan, 12V
EA-160	UV-A Tube, 6 Watts	10.25 in (26.04 cm)	3.25 in (8.26 cm)	2.5 in (6.35 cm)	2.7 lb (1.22 kg)	USA, Australian, United Kingdom, China, Europe, Japan, 12V