

SPECTROLINE®

# E Series

Handheld UV Lamp

EA-140 + EA-160



Forensics



Lab &  
LifeSciences



Mineral  
Inspection

## TOGGLE SWITCH

Built in toggle switch for easy On/Off operation

## ERGONOMIC CARRYING HANDLE

Comfortable grip while operating the lamp.

## ANODIZED ALUMINUM HOUSING

Lightweight and rugged for the toughest work environments.

## BUILT IN REFLECTORS

Corrosion resistant, aluminum reflectors intensify the high UV output.





**EA-140 E SERIES HANDHELD UV LAMP**

**EA-160 E SERIES HANDHELD UV LAMP**



	<b>EA-140</b>	<b>EA-160</b>
<b>Light Source:</b>	UV-A Tube	UV-A Tube
<b>Lamp Style:</b>	Handheld Lamp	Handheld Lamp
<b>Power:</b>	4 Watts	6 Watts
<b>Length:</b>	7 in (17.8 cm)	10.25 in (26.03 cm)
<b>Width:</b>	3.25 in (8.25 cm)	3.25 in (8.25 cm)
<b>Height:</b>	2.5 in (6.35 cm)	2.5in (6.35 cm)
<b>Weight:</b>	2.45 lb (1.1 kg)	2.7 lb (1.22 kg)
<b>Plug Options:</b>	USA, Australian (/FA), United Kingdom (/FBE), China (/FC), European (/FE), Japan(/J), 12 Volt (/12)	



DISTRIBUTED BY:

[www.Spectro-UV.com](http://www.Spectro-UV.com)  
 4 Dubon Ct., Farmingdale, NY 11735  
 866-230-7305