

SPECTROLINE®

E Series

Handheld UV Lamp

ENF-240C + ENF-260C + ENF-280C



Forensics



Lab &
Life Sciences



Mineral
Inspection

ERGONOMIC CARRYING HANDLE
Comfortable grip while operating the lamp.

TOGGLE SWITCH
Built in toggle switch for easy On/Off operation



ANODIZED ALUMINUM HOUSING
Lightweight and rugged for the toughest work environments.

LONG LIFE FILTER GLASS
Our unique filter glass ages at a slower rate than conventional filters and resists solarization extending the lifetime of the lamp.





ENF-240C E SERIES HANDHELD UV LAMP

ENF-260C E SERIES HANDHELD UV LAMP



ENF-280C E SERIES HANDHELD UV LAMP

	ENF-240C	ENF-260C	ENF-280C
Light Source:	4 Watt 365nm Tube / 4 Watt 254nm Tube with filter assembly	4 Watt 365nm Tube / 6 Watt 254nm Tube with filter assembly	8 Watt 365nm Tube / 8 Watt 254nm Tube with filter assembly
Lamp Style:	Handheld Lamp	Handheld Lamp	Handheld Lamp
Length:	8 in (20.32 cm)	10.25 in (26.04 cm)	13.25 in (33.65 cm)
Width:	3.25 in (8.26 cm)	3.25in (8.26cm)	3.25 in (8.26cm)
Height:	2.5 in (6.35 cm)	2.5in (6.35 cm)	2.5 in (6.35cm)
Weight:	2.82 lb (1.28 kg)	3.9 lb (1.77 kg)	1.6 lb (.72 kg)
Plug Options:	USA, Australian (/FA), United Kingdom (/FBE), China (/FC), European (/FE), Japan(/J), 12 Volt (/12H)		