

# Q Series

Blacklight Woods Exam Lamp

Q-12 + Q-22 + Q-22B + QZ-22 + QZ-22B



Forensics



Lab &  
Life Sciences



Mineral  
Inspection

### MAGNIFIER LENS

Capable of 3x magnification at 8" focal length.

### INTEGRATED ALUMINUM REFLECTOR


Provides for maximum irradiance.

### PHENOLIC HANDLE

Rugged handle is attached to the lamp housing by a built in friction lock swivel joint for comfort and flexibility.

### POLYCARBONATE HOUSING

Made from impact resistant polycarbonate housing and ultra-light weight ballast design for fatigue free use.





### Q-12

Typical Peak UV Intensity at 6" (15cm): 700  $\mu$ W/cm<sup>2</sup>

Two 4 Watt Integrally Filtered UV-A Tubes

Net Weight: 2 Lbs (0.9 Kg)



Q-12

### Q-22

Typical Peak UV Intensity at 6" (15cm): 700  $\mu$ W/cm<sup>2</sup>

Two 4 Watt Integrally Filtered UV-A Tubes

Two 4 Watt White Light Tubes

Net Weight: 2 1/4 Lbs (1.0 Kg)



Q-22

### Q-22B

Typical Peak UV Intensity at 6" (15cm): 1400  $\mu$ W/cm<sup>2</sup>

Two 4 Watt Integrally Filtered UV-A Tubes

Net Weight: 2 3/4 Lbs (1.2 Kg)



QZ-22B

### QZ-22

Typical Peak UV Intensity at 6" (15cm): 700  $\mu$ W/cm<sup>2</sup>

Two 4 Watt Integrally Filtered UV-A Tubes

Two 4 Watt White Light Tubes

Net Weight: 2 Lbs (0.9Kg)

### QZ-22B

Typical Peak UV Intensity at 6" (15cm): 1400  $\mu$ W/cm<sup>2</sup>

Four 4 Watt Integrally Filtered UV-A Tubes

Net Weight: 2 Lbs (0.9Kg)



[www.Spectro-UV.com](http://www.Spectro-UV.com)

4 Dubon Ct., Farmingdale, NY 11735

866-230-7305

DISTRIBUTED BY: